

Product Information	
Catalog No.	FC13041
Product name	Anti-Human SLC5A8 Antibody (SAA0314), APC
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	VHH-mFC
Clone ID	SAA0314
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Apical iodide transporter, Sodium-coupled monocarboxylate transporter 1, SMCT1, AIT, Sodium iodide-related cotransporter, SMCT, Solute carrier family 5 member 8, Electrogenic sodium monocarboxylate cotransporter, SLC5A8
Accession	Q8N695
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.