

Product Information	
Catalog No.	FC13053
Product name	Anti-Mouse ENTPD2/CD39L1 Antibody (SAA0770), APC
Applications	ELISA, FCM
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0770
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Ecto-ATP diphosphohydrolase 2, Ecto-ATPDase 2, ENTPD2, CD39 antigen-like 1, Ectonucleoside triphosphate diphosphohydrolase 2, NTPDase 2, Ecto-ATPase 2, CD39L1
Accession	O55026
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.