

Product Information	
Catalog No.	FC13122
Product name	Anti-Human CXCL13/BCA-1/BLC Antibody (SAA0468), PerCP
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, lambda
Clone ID	SAA0468
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	SCYB13, BLC, CXCL13, B cell-attracting chemokine 1, B lymphocyte chemoattractant, Small-inducible cytokine B13, CXC chemokine BLC, Angie, C-X-C motif chemokine 13, BCA-1, BCA1
Accession	O43927
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.