

Product Information	
Catalog No.	FC13229
Product name	Anti-Human CD74 Antibody (SP7219), PerCP
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	SP7219
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	HLA-DR antigens-associated invariant chain, HLA class II histocompatibility antigen gamma chain, CD74, CLIP, Ia antigen-associated invariant chain, DHLAG, Ii
Accession	P04233
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.