

Product Information	
Catalog No.	FC13238
Product name	Anti-Human CYBB/NOX2 Antibody (SAA0847), PerCP
Applications	FCM, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0847
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Cytochrome b558 subunit beta, Cytochrome b(558) subunit beta, Heme-binding membrane glycoprotein gp91phox, Neutrophil cytochrome b 91 kDa polypeptide, gp91-phox, NOX2, Cytochrome b-245 heavy chain, p22 phagocyte B-cytochrome, Superoxide-generating NADPH oxidase heavy chain subunit, gp91-1, NADPH oxidase 2, CYBB, CGD91-phox
Accession	P04839
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.