

Product Information	
Catalog No.	FC13314
Product name	Anti-Human CD41/ITGA2B & CD61/ITGB3 Antibody (OPG2), PerCP
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	OPG2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CD41, ITGAB, GPIIb, ITGA2B, Platelet membrane glycoprotein IIb, GP2B, Integrin alpha-IIb, GPalpha IIb, GPIIIa, Integrin beta-3, CD61, ITGB3, Platelet membrane glycoprotein IIIa, GP3A
Accession	P08514 & P05106
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.