

Product Information	
Catalog No.	FC13318
Product name	Anti-Human CD14 Antibody (F1031), PerCP
Applications	ELISA, FCM, WB
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	F1031
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Myeloid cell-specific leucine-rich glycoprotein, CD14, Monocyte differentiation antigen CD14, RHC33702
Accession	P08571
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.