

Product Information	
Catalog No.	FC13323
Product name	Anti-Mouse CD45/PTPRC Antibody (SAA0570), PerCP
Applications	FCM
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	VHH, mFc (IgG2a)
Clone ID	SAA0570
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	T200, Leukocyte common antigen, CD45, L-CA, Receptor-type tyrosine-protein phosphatase C, PTPRC
Accession	P06800
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.