

Product Information	
Catalog No.	FC13327
Product name	Anti-Human ADRB1 Antibody (6B9), PerCP
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	VHH-hFc
Clone ID	6B9
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Beta-1 adrenoreceptor, ADRB1, Beta-1 adrenergic receptor, ADRB1R, B1AR, Beta-1 adrenoceptor
Accession	P08588
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.