

Product Information	
Catalog No.	FC13328
Product name	Anti-Human CD16 Antibody (3G8), PerCP
Applications	ELISA, FCM, IF
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	3G8
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CD16, CD16a, CD16A, FCGR3A, Fc-gamma RIII-alpha, FcRIIIa, IgG Fc receptor III-2, Low affinity immunoglobulin gamma Fc region receptor III-A, IGFR3, CD16a antigen, Fc-gamma RIII, FCGR3, Fc-gamma RIIIa, FCG3, FcRIII, FcR-10, Low affinity immunoglobulin gamma Fc region receptor III-B, Fc-gamma RIII-beta, CD16-I, Fc-gamma RIII, Fc-gamma RIIIb, FcRIII, FcRIIIb, FcR-10, IgG Fc receptor III-1, CD16b, FCGR3B, CD16B, FCG3, FCGR3, IGFR3
Accession	P08637 & O75015
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.