

Product Information	
Catalog No.	FC13413
Product name	Anti-Human CD121a/IL1R1 Antibody (15C4), PerCP
Applications	ELISA, FCM, IF, WB
Species reactivity	Human, Monkey
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	15C4
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL-1RT1, IL-1R-alpha, mL-1R1, CD121a, IL-1R-1, sIL-1RI, IL1RT1, IL1R1, sIL-1R1, IL1R, CD121 antigen-like family member A, Interleukin-1 receptor type I, Interleukin-1 receptor alpha, mL-1RI, IL1RA, IL-1RT-1, Interleukin-1 receptor type 1, p80
Accession	P14778
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.