

Product Information	
Catalog No.	FC13422
Product name	Anti-Human CD119/IFNGR1 Antibody (SAA0044), PerCP
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0044
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IFN-gamma-R1, Interferon gamma receptor 1, CDw119, Interferon gamma receptor alpha-chain, IFN-gamma receptor 1, IFNGR1, IFN-gamma-R-alpha, CD119
Accession	P15260
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.