

Product Information	
Catalog No.	FC13440
Product name	Anti-Human CD227/MUC1 Antibody (SAA0075), PerCP
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0075
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	PEMT, CA 15-3, H23AG, Peanut-reactive urinary mucin, PEM, EMA, Cancer antigen 15-3, MUC-1, Episialin, MUC1-CT, MUC1-beta, Tumor-associated mucin, CD227, Polymorphic epithelial mucin, KL-6, MUC1-alpha, MUC1, PUM, Krebs von den Lungen-6, MUC1-NT, Tumor-associated epithelial membrane antigen, Mucin-1, Breast carcinoma-associated antigen DF3, Carcinoma-associated mucin
Accession	P15941
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.