

Product Information	
Catalog No.	FC13449
Product name	Anti-Mouse CD140a/PDGFR α Antibody (2.451.1.1b#), PerCP
Applications	ELISA, FCM
Species reactivity	Human, Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	2.451.1.1b#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Platelet-derived growth factor receptor alpha, PDGF-R-alpha, PDGFR-alpha, Alpha platelet-derived growth factor receptor, Alpha-type platelet-derived growth factor receptor, CD140 antigen-like family member A, Platelet-derived growth factor alpha receptor, CD140a, Pdgfra
Accession	P26618
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 μ g per million cells in 100 μ L volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.