

Product Information	
Catalog No.	FC13472
Product name	Anti-Human CD127/IL7R Antibody (MD-707), PerCP
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, Kappa
Clone ID	MD-707
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL7R, Interleukin-7 receptor subunit alpha, IL-7R-alpha, IL-7 receptor subunit alpha, CD127, CDw127, IL-7R subunit alpha, IL-7RA
Accession	P16871
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.