

Product Information	
Catalog No.	FC13499
Product name	Anti-Human Kappa Antibody (HP6156), PerCP
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	HP6156
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Ig kappa chain C region EU, Ig kappa chain C region, Ig kappa chain C region ROY, Immunoglobulin kappa constant, IGKC, Ig kappa chain C region TI, Ig kappa chain C region OU, Ig kappa chain C region CUM, Ig kappa chain C region AG
Accession	P01834
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.