

Product Information	
Catalog No.	FC13506
Product name	Anti-Human CD120b/TNFRSF1B/TNFR2 Antibody (E4#), PerCP
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	E4#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	TNF-R2, TNFR-II, p80 TNF-alpha receptor, Etanercept, TNFBR, TBP-2, TBPII, Tumor necrosis factor receptor superfamily member 1B, CD120b, p75, TNFRSF1B, Tumor necrosis factor receptor 2, TNF-RII, TNFR2, Tumor necrosis factor receptor type II
Accession	P20333
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.