

Product Information	
Catalog No.	FC13538
Product name	Anti-Human CD334/FGFR4 Antibody (U3-17846J6Ye), PerCP
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	U3-17846J6Ye
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	FGFR4, Fibroblast growth factor receptor 4, TKF, FGFR-4, JTK2, CD334
Accession	P22455
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.