

Product Information	
Catalog No.	FC13596
Product name	Anti-Human A2AR/ADORA2A Antibody (SAA0318), PerCP
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0318
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Adenosine receptor A2a, ADORA2A, ADORA2
Accession	P29274
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.