

Product Information	
Catalog No.	FC13617
Product name	Anti-Human CD191/CCR1 Antibody (5908), PerCP
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	5908
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CMKR1, SCYAR1, Macrophage inflammatory protein 1-alpha receptor, CMKBR1, MIP-1alpha-R, HM145, CC-CKR-1, RANTES-R, LD78 receptor, CCR-1, C-C chemokine receptor type 1, CCR1, CD191, C-C CKR-1
Accession	P32246
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.