

Product Information	
Catalog No.	FC13686
Product name	Anti-Human CD178/FASLG/TNFSF6 Antibody (NOK2/NOK-2), PerCP
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	NOK2/NOK-2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Receptor-binding FasL ectodomain, APL, APT1LG1, FasL ICD, Fas antigen ligand, CD178, sFasL, CD95L, APTL, TNFSF6, Apoptosis antigen ligand, SPA, FASLG, FasL, SPPL2A-processed FasL form, Tumor necrosis factor ligand superfamily member 6, Soluble Fas ligand, Fas ligand, CD95 ligand, CD95-L, FASL
Accession	P48023
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.