

Product Information	
Catalog No.	FC13739
Product name	Anti-Mouse CD83 Antibody (SAA2066), PerCP
Applications	FCM
Species reactivity	Human, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA2066
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CD83, hCD83, CD83 antigen, B-cell activation protein, Cell surface protein HB15
Accession	O88324
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.