

Product Information	
Catalog No.	FC13827
Product name	Anti-Human/Mouse/Rat MOG Antibody (8-18C5), PerCP
Applications	ELISA, FCM, IF, WB
Species reactivity	Human, Mouse, Rat, Bovine
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	8-18C5
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	MOG, Myelin-oligodendrocyte glycoprotein
Accession	Q63345
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.