

Product Information	
Catalog No.	FC13844
Product name	Anti-Human NCR3LG1/B7-H6 Antibody (SAA2025), PerCP
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA2025
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	B7H6, NCR3LG1, B7-H6, Natural cytotoxicity triggering receptor 3 ligand 1, B7 homolog 6
Accession	Q68D85
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.