

Product Information	
Catalog No.	FC13879
Product name	Anti-Human CCL17/TARC Antibody (SAA0460), PerCP
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0460
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Thymus and activation-regulated chemokine, SCYA17, TARC, Small-inducible cytokine A17, C-C motif chemokine 17, CC chemokine TARC, CCL17
Accession	Q92583
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.