

Product Information	
Catalog No.	FC13898
Product name	Anti-Human CD268/TNFRSF13C Antibody (VAY-736), PerCP
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	VAY-736
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	BR3, TNFRSF13C, BAFF receptor, BLyS receptor 3, Tumor necrosis factor receptor superfamily member 13C, CD268, B-cell-activating factor receptor, BAFFR, BAFF-R
Accession	Q96RJ3
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.