

Product Information	
Catalog No.	FC13903
Product name	Anti-Human OSMR/IL-31RB Antibody (SAA0162), PerCP
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, Lambda
Clone ID	SAA0162
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	OSMRB, IL-31R subunit beta, Oncostatin-M-specific receptor subunit beta, Interleukin-31 receptor subunit beta, IL-31RB, IL-31R-beta, OSMR, IL-31 receptor subunit beta
Accession	Q99650
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.