

Product Information	
Catalog No.	FC13961
Product name	Anti-Human CD257/BAFF/TNFSF13B Antibody (SAA0081), PerCP
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, kappa
Clone ID	SAA0081
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	BAFF, ZTNF4, BLYS, Tumor necrosis factor ligand superfamily member 13B, Dendritic cell-derived TNF-like molecule, TALL-1, TNF- and APOL-related leukocyte expressed ligand 1, BLYS, B lymphocyte stimulator, TNFSF13B, CD257, TNFSF20, TALL1, B-cell-activating factor, Neutrokine-alpha
Accession	Q9Y275
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.