

Product Information	
Catalog No.	FC13989
Product name	Anti-Human CD294/PTGDR2 Antibody (SAA0325), PerCP
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0325
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Prostaglandin D2 receptor 2, CD294, CRTH2, Chemoattractant receptor-homologous molecule expressed on Th2 cells, GPR44, G-protein coupled receptor 44, DL1R, PTGDR2
Accession	Q9Y5Y4
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.