

Product Information	
Catalog No.	FC13999
Product name	Anti-Human PDGFD Antibody (SAA0451), PerCP
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0451
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IEGF, PDGFD latent form, PDGF-D, Spinal cord-derived growth factor B, SCDGF-B, Iris-expressed growth factor, PDGFD, Platelet-derived growth factor D, PDGFD receptor-binding form, SCDGFB
Accession	Q9GZP0
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.