

Product Information	
Catalog No.	FC14011
Product name	Anti-Mouse TCR $\alpha\beta$ /alpha-beta TCR Antibody (Desire-1), PerCP
Applications	FCM, IF
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	Desire-1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	T cell receptor, T-cell receptor, Ti
Accession	P01852 & P01849
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 μ g per million cells in 100 μ L volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.