

Product Information	
Catalog No.	FC14036
Product name	Anti-Human CD8 Antibody (63H6), PerCP
Applications	FCM
Species reactivity	Human
Clonality	Monoclonal
Isotype	DARPin
Clone ID	63H6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CD8, T-cell surface glycoprotein CD8
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.