

Product Information	
Catalog No.	FC14089
Product name	Anti-Human S1PR3 Antibody (AS55), PerCP
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	AS55
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	EDG3, S1P3, Sphingosine 1-phosphate receptor Edg-3, S1P receptor Edg-3, Sphingosine 1-phosphate receptor 3, S1PR3, Endothelial differentiation G-protein coupled receptor 3, S1P receptor 3
Accession	Q99500
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.