

Product Information	
Catalog No.	FC14092
Product name	Anti-Glycosaminoglycans/GAGs Antibody (SAA2199), PerCP
Applications	ELISA, FCM
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	SAA2199
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Glycosaminoglycans, GAGs
Accession	CAS: 94945-04-7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.