

Product Information	
Catalog No.	FC14114
Product name	Anti-Human CD262/TNFRSF10B/TRAIL-R2 Antibody (SAA1413), PE
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1413
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	TRAILR2, TRAIL-R2, Death receptor 5, CD262, DR5, Tumor necrosis factor receptor superfamily member 10B, ZTNFR9, TRAIL receptor 2, TRICK2, TNFRSF10B, KILLER, TNF-related apoptosis-inducing ligand receptor 2
Accession	O14763
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.