

Product Information	
Catalog No.	FC14127
Product name	Anti-Human CD258/TNFSF14/LIGHT Antibody (SAA0082), PE
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0082
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Herpes virus entry mediator ligand, LIGHT, CD258, HVEM-L, HVEM-L, TNFSF14, Tumor necrosis factor ligand superfamily member 14, Herpesvirus entry mediator ligand
Accession	O43557
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.