

Product Information	
Catalog No.	FC14146
Product name	Anti-Human TNFSF15/TL1A Antibody (320-587), PE
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, lambda
Clone ID	320-587
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	TNFSF15, Vascular endothelial cell growth inhibitor, VEGI, TL1, TNF ligand-related molecule 1, Tumor necrosis factor ligand superfamily member 15, TL1A, TNF-like ligand 1A
Accession	O95150
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.