

Product Information	
Catalog No.	FC14223
Product name	Anti-Mouse CD71/TFRC Antibody (R17217), PE
Applications	FCM
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	R17217
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Transferrin receptor protein 1, TR, TfR, TfR1, Trfr, CD71, Tfrc, R17 217.1.3, R17217, TIB-219
Accession	Q62351
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.