

Product Information	
Catalog No.	FC14258
Product name	Anti-Human IL5 Antibody (SAA0373), PE
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, Kappa
Clone ID	SAA0373
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL5, B-cell differentiation factor I, IL-5, T-cell replacing factor, TRF, Eosinophil differentiation factor, Interleukin-5
Accession	P05113
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.