

Product Information	
Catalog No.	FC14308
Product name	Anti-Human CD271/NGFR Antibody (T-2), PE
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	T-2
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Low affinity neurotrophin receptor p75NTR, TNFRSF16, CD271, Gp80-LNGFR, p75 ICD, Low-affinity nerve growth factor receptor, NGFR, NGF receptor, Tumor necrosis factor receptor superfamily member 16
Accession	P08138
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.