

Product Information	
Catalog No.	FC14484
Product name	Anti-Human CD49b/ITGA2 Antibody (SAA0018), PE
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA0018
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Collagen receptor, ITGA2, VLA-2 subunit alpha, CD49B, CD49 antigen-like family member B, Platelet membrane glycoprotein Ia, GPIa, Integrin alpha-2, CD49b
Accession	P17301
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.