

Product Information	
Catalog No.	FC14520
Product name	Anti-Mouse CD11c/ITGAX Antibody (N418), PE
Applications	FCM, IF
Species reactivity	Mouse
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	N418
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Integrin alpha-X, CD11c, Leukocyte adhesion receptor p150, 95, Leu M5, Leukocyte adhesion glycoprotein p150, 95 alpha chain, CD11 antigen-like family member C, CD11C, ITGAX
Accession	Q9QXH4
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.