

Product Information	
Catalog No.	FC14563
Product name	Anti-Mouse TNC Antibody (81C6), PE
Applications	ELISA, FCM
Species reactivity	Human, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2b, kappa
Clone ID	81C6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	Myotendinous antigen, Neuronectin, GMEM, Tenascin, Tenascin-C, TN-C, TNC, TN, Cytotactin, GP 150-225, Glioma-associated-extracellular matrix antigen, HXB, Hexabrachion, JI
Accession	Q80YX1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.