

Product Information	
Catalog No.	FC14585
Product name	Anti-Human ACVR2A & ACVR2B Antibody (SAA0125), PE
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, lambda
Clone ID	SAA0125
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	ACVR2, ACTRIIA, ACTR-IIA, Activin receptor type IIA, ACVR2A, Activin receptor type-2A, ACTR-IIB, Activin receptor type IIB, Activin receptor type-2B, ACVR2B
Accession	P27037 & Q13705
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.