

Product Information	
Catalog No.	FC14674
Product name	Anti-Mouse CD86/B7-2 Antibody (GL1), PE
Applications	FCM
Species reactivity	Mouse
Host species	Rat
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	GL1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	BU63, B70, CD28LG2, CTLA-4 counter-receptor B7.2, Activation B7-2 antigen, T-lymphocyte activation antigen CD86, CD86, FUN-1
Accession	P42082
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.