

Product Information	
Catalog No.	FC14682
Product name	Anti-Human CD158a/b1/b2/h/j Antibody (SAA0063), PE
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, kappa
Clone ID	SAA0063
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	CD158 antigen-like family member A, KIR2DL1, Killer cell immunoglobulin-like receptor 2DL1, Natural killer-associated transcript 1, NKAT-1, MHC class I NK cell receptor, CD158A, p58 natural killer cell receptor clones CL-42/47.11, CD158a, p58.1 MHC class-I-specific NK receptor, NKAT1, p58 NK receptor CL-42/47.11, p58 NK receptor CL-43, NKAT-6, NKAT6, MHC class I NK cell receptor, CD158 antigen-like family member B1, CD158b1, p58 natural killer cell receptor clone CL-43, Natural killer-associated transcript 6, CD158B1, KIR2DL2, Killer cell immunoglobulin-like receptor 2DL2, P58 natural killer cell receptor clones CL-39/CL-17, CD158i, MHC class I NK cell receptor, CD158I, p58 NK receptor CL-39/CL-17, NKAT8, KKA3, NKAT-8, CD158 antigen-like family member I, Killer cell immunoglobulin-like receptor 2DS4, Natural killer-associated transcript 8, KIR2DS4
Accession	P43626, P43627, P43632
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.