

Product Information	
Catalog No.	FC14741
Product name	Anti-Human CD213a1/IL13RA1 Antibody (SAA0070), PE
Applications	FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA0070
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	IL13R, IL-13R subunit alpha-1, IL13RA1, IL13RA, IL-13RA1, CD213a1, Interleukin-13 receptor subunit alpha-1, IL-13 receptor subunit alpha-1, CT19, IL-13R-alpha-1, Cancer/testis antigen 19
Accession	P78552
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.