

Product Information	
Catalog No.	FC14750
Product name	Anti-Human IL1RL1/ST2 Antibody (SAA1400), PE
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	SAA1400
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	T1, DER4, Interleukin-1 receptor-like 1, ST2, IL1RL1, Protein ST2
Accession	Q01638
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.