

Product Information	
Catalog No.	FC14755
Product name	Anti-Human ROR2/NTRKR2 Antibody (GK-5A1), PE
Applications	ELISA, FCM
Species reactivity	Human, Monkey, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	GK-5A1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	ROR2, Neurotrophic tyrosine kinase, receptor-related 2, NTRKR2, Tyrosine-protein kinase transmembrane receptor ROR2
Accession	Q01974
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.