

Product Information	
Catalog No.	FC14762
Product name	Anti-Human CD202b/TEK Antibody (2339#), PE
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	2339#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.
Concentration	0.5 mg/ml
Target	VMCM, TIE2, Tunica interna endothelial cell kinase, CD202b, Tyrosine-protein kinase receptor TIE-2, TEK, Angiopoietin-1 receptor, VMCM1, Endothelial tyrosine kinase, Tyrosine-protein kinase receptor TEK, p140 TEK, Tyrosine kinase with Ig and EGF homology domains-2, hTIE2
Accession	Q02763
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Store at 4°C for 12 months. Protect from light. Do not freeze.
Note	For flow cytometric staining, the suggested use of this reagent is 0.5 µg per million cells in 100 µL volume. It is recommended that the reagent be titrated for optimal performance for each application. For research use only.